

How to choose the right size

The MIMAR® Suit Recovery size is chosen by the **patient's weight**, not breed. The model, long or short, depends on the area to protect.

Three steps

- 1 Check the weight**
 Use the patient's current weight, from the consultation or clinical record.
- 2 Find the size**
 Locate the weight in the chart. Breeds are only a guide: weight is what counts.
- 3 Choose the model**
 Long suit for abdomen, thorax and back; short suit for neck, thorax and oesophagostomy tubes.

Size chart

Size	Weight	Dog (examples)	Cat (examples)
T1	Up to 4 kg	Chihuahua, Yorkshire	Singapura, Devon Rex
T2	4–10 kg	Poodle, Dachshund	Persian, Ragdoll
T3	10–25 kg	Beagle, Cocker	Maine Coon
T4	25–38 kg	Labrador, Boxer	—
T5	38–48 kg	German Shepherd	—
T6	Over 48 kg COMING SOON	Mastiff, St Bernard	—

Between two sizes?

Try both and choose the one that fits best. Also consider the patient's build when it is atypical: very deep chest (Greyhound, Whippet), very broad (Bulldog) or growing puppies. If in doubt, we can advise you.

How to check the fit

Body, snug: the fabric is very elastic. **At the neck,** with the rigid neck plate in place, leave half a finger of slack: snug but not hurting. If too loose, the animal could pull it off over its head: **always check that it cannot pull its head out through the neck.**



In cats

Most cats wear sizes 1 and 2, but large breeds such as the Maine Coon may need size 3: always choose by weight. Except in special situations, the rigid neck plate is not usually needed: instead, fit the extra Velcro strips included in the bag, from the neck to the back.